

Icancalve Model D



About ICANCLAVE D

ICANCLAVE model D Series is a Class B autoclave in complete compliance with European Standard EN13060. This series includes 8L(normal and quick version), 12L(normal and quick version), 18L, 23L and 29L(super long chamber vesion). They are appropriate for the sterilization of all non-wrapped or wrapped solid, hollow and porous products as represented by the test loads in this standard. ICANCLAVE Model D has been widely used in dentistry, medical, laboratory, beauty, veterinary, pedicure and tattoo industry.




Effiency

Utilisation of an intelligent systems and high performance components (such as instantaneous steam generator and high effecient vacuum pump), which provide a fast sterilization cycle.



Reliability and safety

Full compliance with the European sterilizer standard and certified quality. Ensure a completed safety sterilization result.



Utilizers mechanical and electronic safety mechanisms for pressure & temperature overload protection and safety door locking system.



For any load type

Class B steam sterilizer can be used sterilize any load, such as the instruments of hollow, porous or solid.



Optional internal printer Allows sterilization cycle reports simply printed on thermal paper afeter cycle end.



Data export

Standard equipped with an internal memory able to store 10000+ cycles, and the cycle reports can be automatically saved via the USB driver.

Technical parameters

Order Code	Model	Voltage	Rated Power	Chamber Size	Weight
D208	STE-8-D	220V/110V	1750W	Ф170x320mm	37kg
D212	STE-12-D	220V/110V	1750W	Ф200x360mm	40kg
D218	STE-18-D	220V/110V	1750W	Ф250x350mm	47kg
D223	STE-23-D	220V/110V	1750W	Ф250x450mm	53kg

Sterilization Cycles

Program	Temp. / Pressure	Steri. Time	(Total Time) D208 / D212 / D218 / D223
Solid 134°C	134°C/210Kpa	4min	17-22min 19-24min 21-33min 22-35min
Solid 121°C	121°C/110Kpa	20min	25-35min 25-35min 30-40min 31-40min
Wrapped 134°C	134°C/210Kpa	4min	22-40min 25-40min 25-45min 30-50min
Wrapped 121°C	121°C/110Kpa	20min	38-50min 42-50min 45-60min 48-60min
Textile 134°C	134°C/210Kpa	8min	25-45min 30-45min 31-50min 36-55min
Textile 121°C	121°C/110Kpa	30min	50-60min 52-60min 53-70min 44-70min
Prion 134°C	134°C/210Kpa	18min	35-55min 40-55min 40-60min 46-65min
Liquid 121°C(Option)	121°C/110Kpa	30min	38-50min 40-50min 43-55min 45-60min
B&D Test	134°C/210Kpa	3.5min	10-20min 11-23min 13-25min 16-30min
Helix Test	134°C/210Kpa	3.5min	10-20min 11-23min 13-25min 16-30min
Vacuum Test	/	/	17-20min 18-21min 19-22min 20-25min



*The cycle time is based on the load and with/without pre-heating. *Drying time can be adjusted.



flash water tank



stainless steel chamber



water filter

USB port



draining connector draining connector for used water tank for flash water tank

LCD Screen



Model D Quick series -- D208Q

D208Q is B class autoclaves belong to ICANCLAVE Model D and has same shape with D208, but it with "Quick" program included.

Interm of the sterilizaiton speed, especially the serilization of exposed instruments in the hot state, and has an independent drying program. It is especially suitable for rapid sterilization of dental handpiece and phthalmology and other micro instruments, greatly improving the efficiency of surgery, and making up for shortcomings of small capacity, high price and insufficient drying of cassette serilizers on the market.

Technical parameters

Program	Temp. / Pressure	Steri. Time	Total Time
Quick 134°C	134°C/210Kpa	3.5min	7-12min
Wrapped 134°C	134°C/210Kpa	4min	20-40min
Wrapped 121°C	121°C/110Kpa	20min	38-50min
Prion 134°C	134°C/210Kpa	18min	35-45min
Dry	/	/	1-20min
B&D Test	134°C/210Kpa	3.5min	10-20min
Helix Test	134°C/210Kpa	3.5min	10-20min
Vacuum Test	/	/	17-20min

*The cycle time is based on the load and with/without pre-heating. *Drying time can be adjusted.



QUICK cycle 10 minutes Only

D229---Steam Sterilizer for Extra-long instruments

D229 is specially designed for ultra-long instrument. The depth of the sterilizer is up to 620mm, which effectively solves the problem that ultra-long instruments such as laparoscope can not be sterilized in the ordinary benchtop sterilizers.




Depth of Chamber up to 620mm

				-1-	PO
			5	1	1
1		7		First	1
F D		L	ALL DE LE	11,114,22	1994
Samu	रेत्वमस	HHH	-		Start S
STELLER					
Parts				state and	And I
RUH				appender .	>
1200		- and a state of the			

Tray Dimensions of D229

3Traystotal. Tray #1: 185x570mm Tray #2: 213x570mm Tray #3: 227x570mm

Dimension



Order Code	Chamber Size	W	Н	D	L1	L2	Td	Tw
D208	Ф170x320mm	420mm	370mm	525mm	360mm	330mm	265mm	130mm
D208Q	Ф170x320mm	420mm	370mm	525mm	360mm	330mm	265mm	130mm
D212	Ф200x360mm	420mm	370mm	595mm	400mm	330mm	315mm	160mm
D218	Ф250x350mm	490mm	455mm	600mm	405mm	403mm	285mm	192mm
D223	Ф250x450mm	490mm	455mm	690mm	500mm	403mm	380mm	192mm
D229	Ф250x620mm	490mm	455mm	875mm	690mm	403mm	*	*

*There are three different trays for D229, please refer to the dimension on previous page.



REF ICANCLAVE-21022100

ICANCLAVE reserves the right of modifying characteristics and/or details of the products introduced hereunder, without previous notice. Therefore the products can result lightly modified compared to the pictures.